

## Product datasheet for **HP210066**

### **RASA3 Human qPCR Primer Pair (NM\_007368)**

#### **Product data:**

Product Type:	qPCR Primer Pairs
Gene ID:	22821
Forward Sequence:	TTCCACGCCATCACTGAGTCTG
Reverse Sequence:	GGCGCAAAGAACCCTCAGGAAGA
ACCN:	<u><a href="#">NM_007368</a></u> , <u><a href="#">NM_007368.1</a></u> , <u><a href="#">NM_007368.2</a></u> , <u><a href="#">NM_007368.3</a></u> , <u><a href="#">BC038456</a></u> , <u><a href="#">BC038456.1</a></u> , <u><a href="#">BC047242</a></u> , <u><a href="#">NM_007368.4</a></u>
UniProt ID:	<u><a href="#">Q14644</a></u>
Synonyms:	GAP1IP4BP; GAPIII
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)