

Product datasheet for **HP210075**

AKR7A3 Human qPCR Primer Pair (NM_012067)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22977
Forward Sequence:	CAGGCACTTTGGACTGAGGTTC
Reverse Sequence:	ACATCTCTGCCAGGTATTCCC
ACCN:	NM_012067 , NM_012067.1 , NM_012067.2 , BC042420 , BC042420.1 , BC025709 , BC031562 , BG747063 , BM855386 , BM920315 , NM_012067.3
UniProt ID:	O95154
Synonyms:	AFAR2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)