

Product datasheet for **HP210106**

BACE2 Human qPCR Primer Pair (NM_012105)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	25825
Forward Sequence:	GACGAGAACTCCAGCAGGTCAT
Reverse Sequence:	ACCGATCACCAGCGCATTTGTG
ACCN:	<u>NM_012105</u> , <u>NM_012105.1</u> , <u>NM_012105.2</u> , <u>NM_012105.3</u> , <u>NM_012105.4</u> , <u>BC014453</u> , <u>BC014453.1</u> , <u>BC020582</u> , <u>BQ006998</u> , <u>NM_012105.5</u>
UniProt ID:	<u>Q9Y5Z0</u>
Synonyms:	AEPLC; ALP56; ASP1; ASP21; BAE2; CDA13; CEAP1; DRAP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)