

Product datasheet for **HP210116**

FLASH (CASP8AP2) Human qPCR Primer Pair (NM_012115)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9994
Forward Sequence:	CTCTGAATAGTCCAGTGAGACCT
Reverse Sequence:	CTGACTCTTAGAGGTAGAATGAGC
ACCN:	NM_012115 , NM_012115.1 , NM_012115.2 , NM_012115.3 , BC132828 , BC132830 , BC012574 , BC027850 , BC029386 , BC042577 , BC056685 , BC150274 , BP192859 , BP363682 , NM_012115.4
UniProt ID:	Q9UKL3
Synonyms:	CED-4; FLASH; RIP25
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)