

Product datasheet for **HP210119**

Nocturnin (NOCT) Human qPCR Primer Pair (NM_012118)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	25819
Forward Sequence:	GCAATGGAACATCCTCGCCCAA
Reverse Sequence:	GCCAGGATTTCTTCCAGGATGAG
ACCN:	BC077723 , NM_012118 , NM_012118.1 , NM_012118.2 , NM_012118.3 , BC077723.1 , BC113494 , BC021963 , BC023512 , BC113500 , NM_012118.4
UniProt ID:	Q9UK39
Synonyms:	Ccr4c; CCR4L; CCRN4L; NOC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)