

## Product datasheet for **HP210120**

### CD2AP Human qPCR Primer Pair (NM\_012120)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23607
Forward Sequence:	CCAAAGCCTGAACTGATAGCTGC
Reverse Sequence:	GGACTTGTGGAGCTGCTGGTTT
ACCN:	<u><a href="#">NM_012120</a></u> , <u><a href="#">NM_012120.1</a></u> , <u><a href="#">NM_012120.2</a></u> , <u><a href="#">BC069444</a></u> , <u><a href="#">BC069444.1</a></u> , <u><a href="#">BU679191</a></u> , <u><a href="#">NM_012120.3</a></u>
UniProt ID:	<u><a href="#">Q9Y5K6</a></u>
Synonyms:	CMS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)