

## Product datasheet for **HP210127**

### Claudin 12 (CLDN12) Human qPCR Primer Pair (NM\_012129)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	9069
Forward Sequence:	GAGCAGTGACTGCCTGATGTAC
Reverse Sequence:	GCACATTCCAATCAGGCAGAGC
ACCN:	<u><a href="#">NM_012129</a></u> , <u><a href="#">NM_012129.1</a></u> , <u><a href="#">NM_012129.2</a></u> , <u><a href="#">NM_012129.3</a></u> , <u><a href="#">NM_012129.4</a></u> , <u><a href="#">BC068532</a></u> , <u><a href="#">BC068532.1</a></u> , <u><a href="#">BC036754</a></u> , <u><a href="#">NM_012129.5</a></u>
UniProt ID:	<u><a href="#">P56749</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)