

## Product datasheet for **HP210142**

### DUX4L8 Human qPCR Primer Pair (NM\_012147)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 26583  |
| Forward Sequence: | CGAAGGAGACTCGTTTGGACCC   |
| Reverse Sequence: | CAGCTGGCATGACCTCTCATTC   |
| ACCN:             | <u><a href="#">NM_012147</a></u> , <u><a href="#">NM_012147.1</a></u> , <u><a href="#">NM_012147.2</a></u> , <u><a href="#">NM_012147.3</a></u> , <u><a href="#">NM_012147.4</a></u>   |
| Synonyms:         | DUX2   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)