

## Product datasheet for **HP210154**

### FBXL4 Human qPCR Primer Pair (NM\_012160)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26235
Forward Sequence:	ACAGTCTCGCTGCACTCTTGTC
Reverse Sequence:	GGCTGCAAGACAATTCAAGGCG
ACCN:	<a href="#">NM_012160</a> , <a href="#">NM_012160.1</a> , <a href="#">NM_012160.2</a> , <a href="#">NM_012160.3</a> , <a href="#">NM_012160.4</a> , <a href="#">BC055010</a> , <a href="#">BC055010.1</a> , <a href="#">BC032641</a> , <a href="#">BC039819</a> , <a href="#">BC091484</a> , <a href="#">BQ213738</a>
UniProt ID:	<a href="#">Q9UKA2</a>
Synonyms:	FBL4; FBL5; MTDPS13
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)