

Product datasheet for **HP210197**

OGA Human qPCR Primer Pair (NM_012215)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10724
Forward Sequence:	GCAAGAGTTTGGTGTGCCTCATC
Reverse Sequence:	GTGCTGCAACTAAAGGAGTCCC
ACCN:	<u>NM_012215</u> , <u>NM_012215.1</u> , <u>NM_012215.2</u> , <u>NM_012215.3</u> , <u>BC001343</u> , <u>BC014201</u> , <u>BC032591</u> , <u>BC039583</u> , <u>BC047877</u> , <u>BQ016159</u> , <u>NM_012215.5</u>
UniProt ID:	<u>O60502</u>
Synonyms:	MEA5; MGEA5; NCOAT
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)