

Product datasheet for **HP210201**

NEK1 Human qPCR Primer Pair (NM_012224)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4750
Forward Sequence:	GGACAGTATGAACATTACCATGCC
Reverse Sequence:	CTGGACGAACTCCAGGCAGAAT
ACCN:	NM_012224 , NM_012224.1 , BC015147 , BC037790 , BC068201 , BC114491 , NM_012224.3 , NM_012224.4
UniProt ID:	Q96PY6
Synonyms:	ALS24; NY-REN-55; SRPS2; SRPS2A; SRTD6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)