

Product datasheet for **HP210202**

NUBP2 Human qPCR Primer Pair (NM_012225)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10101
Forward Sequence:	TGCACCAGCGTCTTCTCCAGG
Reverse Sequence:	GAACTCCTGGATGAAGTCGTGG
ACCN:	<u>NM_012225</u> , <u>NM_012225.1</u> , <u>NM_012225.2</u> , <u>NM_012225.3</u> , <u>BC008005</u> , <u>BC002768</u> , <u>BF439443</u> , <u>BI196883</u> , <u>BU189242</u> , <u>BU680530</u> , <u>BX446436</u> , <u>BX648907</u> , <u>NM_012225.4</u>
UniProt ID:	<u>Q9Y5Y2</u>
Synonyms:	CFD1; CIAO6; NBP 2; NUBP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)