

## Product datasheet for **HP210208**

### CAVIN1 Human qPCR Primer Pair (NM\_012232)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	284119
Forward Sequence:	GGCAGATCAAGAAGCTGGAGGT
Reverse Sequence:	CAGCGATTGCTGATGCTCAGTT
ACCN:	<a href="#">NM_012232</a> , <a href="#">NM_012232.1</a> , <a href="#">NM_012232.2</a> , <a href="#">NM_012232.3</a> , <a href="#">NM_012232.4</a> , <a href="#">NM_012232.5</a> , <a href="#">BC066123</a> , <a href="#">BC066123.1</a> , <a href="#">BC004295</a> , <a href="#">BC008849</a> , <a href="#">BC073759</a> , <a href="#">NM_012232.6</a>
UniProt ID:	<a href="#">Q6NZI2</a>
Synonyms:	CAVIN; cavin-1; CGL4; FKSG13; PTRF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)