

Product datasheet for **HP210209**

RAB3GAP1 Human qPCR Primer Pair (NM_012233)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22930
Forward Sequence:	GAGTGCCTGTCTGCTCTCAGAT
Reverse Sequence:	CTCCTTTCAGCACCATTTGCC
ACCN:	NM_012233 , NM_012233.1 , NM_012233.2 , BC020303 , BC022977 , BC030992 , BC036732 , BC071602 , BC146809 , BP872264 , NM_012233.3
UniProt ID:	Q15042
Synonyms:	P130; RAB3GAP; RAB3GAP130; WARBM1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)