

Product datasheet for **HP210228**

HEY2 Human qPCR Primer Pair (NM_012259)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23493
Forward Sequence:	TGAGAAGACTTGTGCCAACTGCT
Reverse Sequence:	CCCTGTTGCCTGAAGCATCTTC
ACCN:	BC007707 , NM_012259 , NM_012259.1 , NM_012259.2 , BC007707.2 , BI549391 , BT020067 , BT020068 , NM_012259.3
UniProt ID:	Q9UBP5
Synonyms:	bHLHb32; CHF1; GRIDLOCK; GRL; HERP1; HESR2; HRT2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)