

Product datasheet for **HP210234**

DNAJB5 Human qPCR Primer Pair (NM_012266)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	25822
Forward Sequence:	GTGCTCTACAGTGCCTGATCA
Reverse Sequence:	TCTTCACGGTGCCTGGCTTGAT
ACCN:	NM_012266 , NM_012266.1 , NM_012266.2 , NM_012266.3 , NM_012266.4 , NM_012266.5 , BI915181 , BU845478 , BX419166 , NM_012266.6
UniProt ID:	O75953
Synonyms:	Hsc40
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)