

Product datasheet for **HP210236**

HYAL4 Human qPCR Primer Pair (NM_012269)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23553
Forward Sequence:	CACCATAGGAGAAAGTGCTGCC
Reverse Sequence:	TCAGCAGCTCTGGTCACATTGG
ACCN:	NM_012269 , NM_012269.1 , NM_012269.2 , BC104788 , BC104788.1 , BC104790 , BC143654 , BF085228 , BG013757 , NM_012269.3
UniProt ID:	Q2M3T9
Synonyms:	CSHY; HYAL-4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)