

Product datasheet for **HP210239**

IL36RN Human qPCR Primer Pair (NM_012275)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26525
Forward Sequence:	TGGAAGCCAGTGCCTGTCATGT
Reverse Sequence:	AGCCGACTCGAAGCTGGAGGT
ACCN:	<u>NM_012275, NM_012275.1, NM_012275.2, BC024747, BC000945, NM_012275.3</u>
UniProt ID:	<u>Q9UBH0</u>
Synonyms:	FIL1; FIL1(DELTA); FIL1D; IL-36Ra; IL1F5; IL1HY1; IL1L1; IL1RP3; IL36RA; PSORP; PSORS14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)