

Product datasheet for **HP210243**

FTSJ1 Human qPCR Primer Pair (NM_012280)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	24140
Forward Sequence:	TGACACGGGCAGTTGACCTGTG
Reverse Sequence:	CACACCTGGTAGTGGAGCCATA
ACCN:	<u>NM_012280</u> , <u>NM_012280.1</u> , <u>NM_012280.2</u> , <u>NM_012280.3</u> , <u>BC023584</u> , <u>BC023584.1</u> , <u>BC021047</u> , <u>NM_012280.4</u>
UniProt ID:	<u>Q9UET6</u>
Synonyms:	CDLIV; JM23; MRX9; MRX44; SPB1; TRMT7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)