

Product datasheet for **HP210255**

HA1 (ARHGAP45) Human qPCR Primer Pair (NM_012292)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23526
Forward Sequence:	GCTTCAACGTGAGTGATGTGGC
Reverse Sequence:	GCCATGACTTGTGAACCTGGTG
ACCN:	NM_012292 , NM_012292.1 , NM_012292.2 , NM_012292.3 , NM_012292.4 , BC065223 , BC065223.1 , BC035564 , BC048129 , BE279735 , BU838153 , BX119705
UniProt ID:	Q92619
Synonyms:	HA-1; HLA-HA1; HMHA1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)