

Product datasheet for **HP210269**

FBXL11 (KDM2A) Human qPCR Primer Pair (NM_012308)

Product data:

| | |
|--------------------------|---|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 22992 |
| Forward Sequence: | CAAGGAGAGTGTGGTGGTTTGCC |
| Reverse Sequence: | ACCTCTCCACAGAGGGAACATG |
| ACCN: | <u>NM_012308</u> , <u>NM_012308.1</u> , <u>NM_012308.2</u> , <u>BC001203</u> , <u>BC027467</u> , <u>BC047371</u> , <u>BC047486</u> , <u>BC056905</u> , <u>BC064360</u> , <u>BE675843</u> , <u>NM_012308.3</u> |
| UniProt ID: | <u>Q9Y2K7</u> |
| Synonyms: | CXXC8; FBL7; FBL11; FBXL11; JHDM1A; LILINA |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)