

## Product datasheet for **HP210281**

### LSM4 Human qPCR Primer Pair (NM\_012321)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	25804
Forward Sequence:	GGAGACGTACAATGGACACCTG
Reverse Sequence:	ATGCGCAGGTACTTGATGGTGC
ACCN:	<a href="#">NM_012321</a> , <a href="#">NM_012321.2</a> , <a href="#">NM_012321.3</a> , <a href="#">NM_012321.4</a> , <a href="#">BC000387</a> , <a href="#">BC003652</a> , <a href="#">BC022198</a> , <a href="#">BC023665</a> , <a href="#">BF432099</a> , <a href="#">BQ069184</a> , <a href="#">BU539495</a> , <a href="#">BX096089</a> , <a href="#">NM_012321.5</a>
UniProt ID:	<a href="#">Q9Y4Z0</a>
Synonyms:	GRP; YER112W
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)