

Product datasheet for **HP210293**

MYO10 Human qPCR Primer Pair (NM_012334)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 4651 |
| Forward Sequence: | CACTCTGCCGTATTTCCACAGC |
| Reverse Sequence: | TTTGTGGAGCCAGCCTTGCTTG |
| ACCN: | NM_012334 , NM_012334.1 , NM_012334.2 , BC014956 , BC032669 , BC041694 , BC050682 , BC071965 , BC094838 , BC108736 , BC137168 , BC150285 |
| UniProt ID: | Q9HD67 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)