

## Product datasheet for **HP210326**

### Pallidin (BLOC1S6) Human qPCR Primer Pair (NM\_012388)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26258
Forward Sequence:	AGTGGAGCAACTGGCAGAAGGA
Reverse Sequence:	ATCTCTTGTTCCAGTGTGTCTAAC
ACCN:	<a href="#">NM_012388</a> , <a href="#">NM_012388.1</a> , <a href="#">NM_012388.2</a> , <a href="#">NM_012388.3</a> , <a href="#">BC004819</a> , <a href="#">BC004819.1</a> , <a href="#">BC026289</a> , <a href="#">BP423071</a> , <a href="#">NM_012388.4</a>
UniProt ID:	<a href="#">Q9UL45</a>
Synonyms:	BLOS6; HPS9; PA; PALLID; PLDN
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)