

Product datasheet for **HP210337**

PLA2G2D Human qPCR Primer Pair (NM_012400)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26279
Forward Sequence:	GCTGCTATGACCACCTGAAGAC
Reverse Sequence:	TGCTCACACCAGCTTCCCTGT
ACCN:	<u>NM_012400</u> , <u>NM_012400.1</u> , <u>NM_012400.2</u> , <u>NM_012400.3</u> , <u>BC025706</u> , <u>BC025706.1</u> , <u>NM_012400.4</u>
UniProt ID:	<u>Q9UNK4</u>
Synonyms:	PLA2IID; sPLA2-IID; sPLA2S; SPLASH
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)