

## Product datasheet for **HP210338**

### Plexin B2 (PLXNB2) Human qPCR Primer Pair (NM\_012401)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23654
Forward Sequence:	CTTCAGCCTGATCCAGAGGTTTG
Reverse Sequence:	GTGGAACACGTAGTCTGTACCC
ACCN:	<u><a href="#">NM_012401</a></u> , <u><a href="#">NM_012401.1</a></u> , <u><a href="#">NM_012401.2</a></u> , <u><a href="#">NM_012401.3</a></u> , <u><a href="#">BC000221</a></u> , <u><a href="#">BC004542</a></u> , <u><a href="#">BC110366</a></u> , <u><a href="#">BC127770</a></u> , <u><a href="#">BC148225</a></u> , <u><a href="#">BC148280</a></u> , <u><a href="#">BC172440</a></u> , <u><a href="#">BT006887</a></u> , <u><a href="#">NM_012401.4</a></u>
UniProt ID:	<u><a href="#">O15031</a></u>
Synonyms:	dj402G11.3; MM1; Nbla00445; PLEXB2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)