

## Product datasheet for **HP210355**

### RLF Human qPCR Primer Pair (NM\_012421)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6018
Forward Sequence:	ACCTTACCCAGCAGCTCCAAAC
Reverse Sequence:	ACACCAAAGTACGACAGCGCT
ACCN:	<u>NM_012421, NM_012421.1, NM_012421.2, NM_012421.3, BC113666, BM542854, NM_012421.4</u>
UniProt ID:	<u>Q13129</u>
Synonyms:	ZN-15L; ZNF292L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)