

Product datasheet for **HP210359**

SF3B3 Human qPCR Primer Pair (NM_012426)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23450
Forward Sequence:	GTGTGCAGTGAACCAGCGACAA
Reverse Sequence:	CTGCTGACATCTCCTTCCGTTC
ACCN:	<u>NM_012426</u> , <u>NM_012426.1</u> , <u>NM_012426.2</u> , <u>NM_012426.3</u> , <u>NM_012426.4</u> , <u>BC068974</u> , <u>BC000463</u> , <u>BC003146</u> , <u>BC009780</u> , <u>BI465052</u>
UniProt ID:	<u>Q15393</u>
Synonyms:	RSE1; SAP130; SF3b130; STAF130
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)