

Product datasheet for **HP210361**

SEC14 like protein 2 (SEC14L2) Human qPCR Primer Pair (NM_012429)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23541
Forward Sequence:	CCTGAAGACCAAGATGGGAGAG
Reverse Sequence:	GCTGTAGGTGTTGTCAAACCGC
ACCN:	NM_012429 , NM_012429.1 , NM_012429.2 , NM_012429.3 , NM_012429.4 , BC058915 , BC048337 , BM680990 , BU734336 , NM_012429.5
UniProt ID:	O76054
Synonyms:	C22orf6; SPF; TAP; TAP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)