

Product datasheet for **HP210387**

TSPAN16 Human qPCR Primer Pair (NM_012466)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	26526
Forward Sequence:	TGCATCACGGTACTGCTTGGCT
Reverse Sequence:	AAAGAAGGACCACTGTGGCAGC
ACCN:	<u>BC029908, NM_012466, NM_012466.1, NM_012466.3, BC029908.1, BM560276, NM_012466.4</u>
UniProt ID:	<u>Q9UKR8</u>
Synonyms:	TM-8; TM4-B; TM4SF16
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)