

Product datasheet for **HP210393**

DNAAF11 Human qPCR Primer Pair (NM_012472)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23639
Forward Sequence:	CGAGCCAACTCAAGGAAGAGG
Reverse Sequence:	TGCCTGTAGGTGGTCTTTGCTC
ACCN:	NM_012472 , NM_012472.1 , NM_012472.2 , NM_012472.3 , NM_012472.4 , NM_012472.5 , BC047286 , BC047286.1 , BC027589 , BG204280 , BM981867 , BU674822 , NM_012472.6
UniProt ID:	Q86X45
Synonyms:	CILD19; LRRC6; LRTP; tiIB; TSLRP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)