

Product datasheet for **HP210418**

VPS4A Human qPCR Primer Pair (NM_013245)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27183
Forward Sequence:	GCCAGTCAAAGAGAACCAGAGTG
Reverse Sequence:	ACGTCGTTCCACCGTATGTTGG
ACCN:	<u>NM_013245</u> , <u>NM_013245.1</u> , <u>NM_013245.2</u> , <u>BC047932</u> , <u>BC047932.1</u> , <u>BC013021</u> , <u>BC035121</u> , <u>NM_013245.3</u>
UniProt ID:	<u>Q9UN37</u>
Synonyms:	CIMDAG; SKD1; SKD1A; SKD2; VPS4; VPS4-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)