

Product datasheet for **HP210419**

NNT1 (CLCF1) Human qPCR Primer Pair (NM_013246)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23529
Forward Sequence:	CCGCTACCTGGAGCACCAACT
Reverse Sequence:	CACACCTCCAAGTCAACAGTGG
ACCN:	NM_013246 , NM_013246.1 , NM_013246.2 , BC012939 , BC066229 , BC066230 , BC066231 , BM846622 , NM_013246.3
UniProt ID:	Q9UBD9
Synonyms:	BSF-3; BSF3; CISS2; CLC; NNT-1; NNT1; NR6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)