

Product datasheet for **HP210425**

TBK1 Human qPCR Primer Pair (NM_013254)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29110
Forward Sequence:	CAACCTGGAAGCGGCAGAGTTA
Reverse Sequence:	ACCTGGAGATAATCTGCTGTCGA
ACCN:	<u>BC009864</u> , <u>NM_013254</u> , <u>NM_013254.1</u> , <u>NM_013254.2</u> , <u>NM_013254.3</u> , <u>BC034950</u> , <u>BC034950.2</u> , <u>NM_013254.4</u>
UniProt ID:	<u>Q9UHD2</u>
Synonyms:	FTDALS4; IIAE8; NAK; T2K
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)