

Product datasheet for **HP210435**

CTNNA3 Human qPCR Primer Pair (NM_013266)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29119
Forward Sequence:	CACCAGAACCTCAGGCAGCAAC
Reverse Sequence:	GGCGTTTCTCTAGTGATGGTCG
ACCN:	BC022004 , NM_013266 , NM_013266.1 , NM_013266.2 , NM_013266.3 , BC030646 , BC065819 , BP362015 , BX103986 , NM_013266.4
UniProt ID:	Q9UI47
Synonyms:	ARVD13; VR22
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)