

Product datasheet for **HP210436**

GLS2 Human qPCR Primer Pair (NM_013267)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27165
Forward Sequence:	TGAGGCACTGTGCTCGGAAGTT
Reverse Sequence:	TCGAAGAGCTGAGACATCGCCA
ACCN:	NM_013267 , NM_013267.1 , NM_013267.2 , NM_013267.3 , BC025661 , BC036929 , BC048344 , BC059370 , BC119767 , BC119768 , BC166649 , BI916699 , NM_013267.4
UniProt ID:	Q9UI32
Synonyms:	GA; GLS; hLGA; LGA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)