

Product datasheet for **HP210441**

DNA Polymerase lambda (POLL) Human qPCR Primer Pair (NM_013274)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	27343
Forward Sequence:	CCATAAGCCTGTCACCTCGTAC
Reverse Sequence:	GCTCTCACTGATATGGTCCAGC
ACCN:	<u>NM_013274</u> , <u>NM_013274.1</u> , <u>NM_013274.2</u> , <u>NM_013274.3</u> , <u>BC068529</u> , <u>BC068529.1</u> , <u>BC003548</u> , <u>BI906881</u>
UniProt ID:	<u>Q9UGP5</u>
Synonyms:	BETAN; POLKAPPA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)