

Product datasheet for **HP210457**

TRA2A Human qPCR Primer Pair (NM_013293)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29896
Forward Sequence:	AATCGGAGAGCAGGTCAGGATC
Reverse Sequence:	GAGAATGTCTCCTTGACCTCGAC
ACCN:	NM_013293 , NM_013293.1 , NM_013293.2 , NM_013293.3 , NM_013293.4 , BC017094 , BC017094.2 , BC069011 , BG722500 , BP376263 , NM_013293.5
UniProt ID:	Q13595
Synonyms:	AWMS1; HSU53209
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)