

Product datasheet for **HP210468**

SLC30A4 Human qPCR Primer Pair (NM_013309)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7782
Forward Sequence:	AAGCCAGGTTGACCATTGCTGC
Reverse Sequence:	TGATGGCGCTTAGGTCAGTTAAC
ACCN:	NM_013309 , NM_013309.1 , NM_013309.2 , NM_013309.3 , NM_013309.4 , NM_013309.5 , BC026089 , BC026089.2 , BC042520 , BM789500 , NM_013309.6
UniProt ID:	O14863
Synonyms:	znT-4; ZNT4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)