

Product datasheet for **HP210472**

PRRC2B Human qPCR Primer Pair (NM_013318)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84726
Forward Sequence:	TGAGAGCCAACAGAATGGGACG
Reverse Sequence:	TTACTGGAGGCTTCGGGAAGGT
ACCN:	BC002872 , NM_013318 , NM_013318.3 , BC002872.2 , BC002969 , BC012289 , BC012451 , BX647528 , BX647842
UniProt ID:	Q5JSZ5
Synonyms:	BAT2L; BAT2L1; KIAA0515; LQFBS-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)