

Product datasheet for **HP210486**

GMPPA Human qPCR Primer Pair (NM_013335)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29926
Forward Sequence:	GGACAGTGAGAGCCTCTTCAAG
Reverse Sequence:	TCGAGTTCAGGATGAGCACCTC
ACCN:	NM_013335 , NM_013335.1 , NM_013335.2 , NM_013335.3 , BC007456 , BC007456.2 , BC014266 , BP284488 , BP366441 , NM_013335.4
UniProt ID:	Q96IJ6
Synonyms:	AAMR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)