

Product datasheet for **HP210489**

ALG5 Human qPCR Primer Pair (NM_013338)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29880
Forward Sequence:	GAGAAGCAGCTTCACGGACGTT
Reverse Sequence:	GTCCAGTTGACAGCAATTCTGC
ACCN:	NM_013338 , NM_013338.1 , NM_013338.2 , NM_013338.3 , NM_013338.4 , BC012531 , BC012531.1 , BF668393 , BG482604 , BP246765 , BU739112 , NM_013338.5
UniProt ID:	Q9Y673
Synonyms:	bA421P11.2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)