

## Product datasheet for **HP210498**

### **NENF Human qPCR Primer Pair (NM\_013349)**

#### **Product data:**

|                          |                                                                                                                                                                                                                                                                                                                                              |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b>     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                            |
| <b>Gene ID:</b>          | 29937                                                                                                                                                                                                                                                                                                                                        |
| <b>Forward Sequence:</b> | TCAGCCCATCTACTTGGCAGTG                                                                                                                                                                                                                                                                                                                       |
| <b>Reverse Sequence:</b> | AGTGTCATGGGTGAGGTCTGCA                                                                                                                                                                                                                                                                                                                       |
| <b>ACCN:</b>             | <u><a href="#">NM_013349</a></u> , <u><a href="#">NM_013349.1</a></u> , <u><a href="#">NM_013349.2</a></u> , <u><a href="#">NM_013349.3</a></u> , <u><a href="#">NM_013349.4</a></u> , <u><a href="#">BC008823</a></u> ,<br><u><a href="#">BC008823.2</a></u> , <u><a href="#">NM_013349.5</a></u>                                           |
| <b>UniProt ID:</b>       | <u><a href="#">Q9UMX5</a></u>                                                                                                                                                                                                                                                                                                                |
| <b>Synonyms:</b>         | CIR2; SCIRP10; SPUF                                                                                                                                                                                                                                                                                                                          |
| <b>Components:</b>       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.                                                                                                                                               |
| <b>Quality Control:</b>  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| <b>Storage:</b>          | Store at -20°C.                                                                                                                                                                                                                                                                                                                              |
| <b>Stability:</b>        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                 |



[View online »](#)