

Product datasheet for **HP210519**

ALIX (PDCD6IP) Human qPCR Primer Pair (NM_013374)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10015
Forward Sequence:	GCTCAGATGAGAGAAGCCACCA
Reverse Sequence:	AGTCTGGATGCCTCCCTGTTCA
ACCN:	<u>NM_013374</u> , <u>NM_013374.2</u> , <u>NM_013374.3</u> , <u>NM_013374.4</u> , <u>NM_013374.5</u> , <u>BC020066</u> , <u>BC020066.1</u> , <u>BC068454</u> , <u>BI492166</u> , <u>BI562417</u> , <u>BM669193</u> , <u>BQ024882</u> , <u>BQ895877</u> , <u>BT007367</u> , <u>BU165197</u> , <u>BU194828</u> , <u>NM_013374.6</u>
UniProt ID:	<u>Q8WUM4</u>
Synonyms:	AIP1; ALIX; DRIP4; HP95
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)