

Product datasheet for **HP210524**

DPP2 (DPP7) Human qPCR Primer Pair (NM_013379)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29952
Forward Sequence:	TTCGCCAGCAACAATGTGACCG
Reverse Sequence:	ATGTTGCTGGCGGCTCTGAGAT
ACCN:	NM_013379 , NM_013379.1 , NM_013379.2 , BC011907 , BC011907.2 , BC016961 , BF342300 , BX400716 , NM_013379.3
UniProt ID:	Q9UHL4
Synonyms:	DPP2; DPPII; QPP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)