

Product datasheet for **HP210526**

POMT2 Human qPCR Primer Pair (NM_013382)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29954
Forward Sequence:	CCGAGTCTATCTGCTTGGCAAC
Reverse Sequence:	CCTCTCTGCATGGCTACAGCAA
ACCN:	NM_013382 , NM_013382.1 , NM_013382.2 , NM_013382.3 , NM_013382.4 , NM_013382.5 , BC031651 , BC031651.1 , BE797202 , BM127609 , BM724049 , BX248027 , NM_013382.7
UniProt ID:	Q9UKY4
Synonyms:	LGMD2N; LGMDR14; MDDGA2; MDDGB2; MDDGC2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)