

Product datasheet for **HP210548**

Ubiquilin (UBQLN1) Human qPCR Primer Pair (NM_013438)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	29979
Forward Sequence:	AGTGGCACTTCTGGGCAGAGTA
Reverse Sequence:	GCAACAAGCTCTGCATTCCTGG
ACCN:	NM_013438 , NM_013438.1 , NM_013438.2 , NM_013438.3 , NM_013438.4 , BC010066 , BC010066.2 , BC039294 , BC039294.1 , NM_013438.5
UniProt ID:	Q9UMX0
Synonyms:	DA41; DSK2; PLIC-1; UBQN; XDRP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)