

Product datasheet for **HP210555**

MKRN1 Human qPCR Primer Pair (NM_013446)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23608
Forward Sequence:	CCAATGGATGCTGCCAGAGAT
Reverse Sequence:	GGTTGGCTTTCTCATAGACCACC
ACCN:	NM_013446 , NM_013446.1 , NM_013446.2 , NM_013446.3 , BC025955 , BC025955.1 , BC012398 , BC037400 , BC051821 , BC064838 , BC127216
UniProt ID:	Q9UHC7
Synonyms:	RNF61
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)