

## Product datasheet for **HP210573**

### DHX38 Human qPCR Primer Pair (NM\_014003)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 9785   |
| Forward Sequence: | GACCTGGATCACTACAGTGCCA   |
| Reverse Sequence: | GTGGCTGATGTGACGATGAGCT   |
| ACCN:             | <a href="#">NM_014003</a> , <a href="#">NM_014003.1</a> , <a href="#">NM_014003.2</a> , <a href="#">NM_014003.3</a> , <a href="#">BC004235</a> , <a href="#">BC004235.2</a> , <a href="#">BC008340</a> , <a href="#">BG719187</a> , <a href="#">BM473438</a> , <a href="#">BU733805</a> , <a href="#">NM_014003.4</a>                        |
| UniProt ID:       | <a href="#">Q92620</a>   |
| Synonyms:         | DDX38; PRP16; PRPF16; RP84   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)