

Product datasheet for **HP210597**

C11orf54 Human qPCR Primer Pair (NM_014039)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	28970
Forward Sequence:	GCTGATGTCCAGGTCTCTGTAG
Reverse Sequence:	AAGGCACACCTCCAATTCTGC
ACCN:	<u>NM_014039</u> , <u>NM_014039.1</u> , <u>NM_014039.2</u> , <u>NM_014039.3</u> , <u>BC012298</u> , <u>BC012298.1</u> , <u>BC007110</u> , <u>BC046477</u> , <u>BI597310</u> , <u>BM551310</u> , <u>NM_014039.4</u>
UniProt ID:	<u>Q9H0W9</u>
Synonyms:	PTD012; PTOD012
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)